


# GIVE YOUR PROTEIN MICROARRAY A BOOST




XanTec's HC Protein microarray slides were designed for superior performance as protein chip substrates. Their preactivated surfaces are coated with a highly protein compatible hydrogel with properties similar to bioinert chromatographic separation media. These hydrogel coatings are, therefore, the optimal choice for sensitive and selective protein detection, even in complex sample matrices like serum, cell lysates, environmental samples or fermentation broths.

The HC coating's unique extended brush structure has a thickness of 1 – 2  $\mu\text{m}$  and allows immobilization levels of  $> 50 \text{ ng/mm}^2$  - up until now only achievable with thick film coatings such as nitrocellulose or polyacrylamide pads. However, unlike these state-of-the-art 3D substrates, the revolutionary HC coatings minimize autofluorescence: as the background of HC hydrogel matrix is extremely low, signal/noise ratios of  $> 1000$  can be achieved. As more than 50% of the spotted protein is covalently immobilized, spotting of low protein quantities below 100 amol/spot is possible, which makes arraying of expensive biomaterial economic. The biocompatible character of the HC hydrogel coating stabilizes the immobilized proteins and suppresses nonspecific binding. Blocking is not required, which significantly enhances the biospecificity of the resulting protein arrays.

## Key features of hydrogel coated HC microarray slides:



- High immobilization capacity of  $> 50 \text{ ng / mm}^2$  ensures maximal signal intensity
- Protein compatible hydrogel matrix preserves activity of immobilized ligands
- Minimal autofluorescence and extremely low level of non-specific binding yield superior S/N ratios
- No blocking step required maximizes specificity of array
- NHS mediated covalent immobilization provides stable binding of proteins, peptides, small molecules, oligonucleotides.
- High immobilization efficiency saves precious biomaterial
- Excellent bioinertness allows work with serum and ultraspecific cell separation



## Featured application of HC hydrogel slides: Cell separation and immobilization

Initially developed for enhancing the sensitivity of protein microarrays, the brush-like structured HC hydrogel coatings are also an excellent surface for selectively attaching cells to solid supports. The non-modified form is inert, consisting of a highly hydrated, slightly negatively charged polymer matrix that effectively prevents non-specific cell adhesion. When coupled to specific receptor molecules (proteins, peptides, carbohydrates or small molecules), a high density of binding sites is generated which, in contrast to normal cross-linked hydrogel matrices, is sterically accessible for the relatively voluminous cells.

The flexible hydrophilic polymer chains of this 1-2  $\mu\text{m}$  thin layer align themselves smoothly to the cell surface and allow the receptor molecules to attach at multiple sites simultaneously. This results in specifically bound cells which are safely embedded in a cell friendly environment. In contrast to common panning techniques, it is not necessary to carefully optimize the incubation parameters or the reaction time as practically no time dependent non-specific adhesion occurs.



**Left Figure:** MDA breast tumour cells (surface antigen EGFR) in human serum were incubated on a partially antibody derivatized HC hydrogel surface.

**Left:** Anti EGFR IgG functionalized HC hydrogel surface.

**Upper right corner:** Blank HC hydrogel surface (control).

For more information contact us:

**XanTec bioanalytics GmbH**

Merowingerplatz 1a  
 D-40225 Düsseldorf  
 Germany

Phone +49 211 993 647 44  
 Fax +49 211 993 647 46  
 E-mail [info@xantec.com](mailto:info@xantec.com)

[www.xantec.com](http://www.xantec.com)

*The above separation technology based on HC hydrogel coatings is winner of the Prize for Medical Technology 2005 & 2009 awarded by the German Federal Ministry of Education and Research.*